Name	Date
Partner	Mr. Romano
AP Environmental Science	

## PPM and PPB HAGD and MTFBWY

## Introduction:

Sometimes people are concerned with materials in very small amounts or concentrations. An environmental scientist might have to analyze water, air, or soil to determine what substances are present, but the amount of the substance is so small that a special unit is required to express it. Well, a unit does exist; just as "per cent" means "out of a hundred", parts per million (ppm) means out of a million and parts per billion (ppb) means out of a billion.

A concentration of 1 ppm corresponds to one part substance per one million parts of the gas, liquid, or solid medium it is found in. If you have any understanding of metric system units: 1 ppm = 1mg/L (not sure if that really helps ...). But, just because these amounts are numerically small, it does not mean that they are insignificant. For example, fish have specific requirements for dissolved oxygen (DO)in water to be able to survive. The DO level must be in the range of 2 and 9 ppm and can be disturbed by environmental activities. When testing groundwater, benzene levels (a carcinogen found in gasoline), are analyzed. The maximum allowable concentration of benzene in groundwater is 1 ppb.

As technology develops, people are able to use more sensitive equipment and testing techniques to detect tiny concentrations of contaminants. The smallest amount that can be detected with the analytical tools available (in this case, your eyes) is known as the *detection limit*.

## **Procedure to Conduct a Successive Serial Dilution:**

1. Place a Chemplate or	top of a clean white piece of pa	per. This will give you a clearer	view and better contrast.	
•	ood coloring vertically in order to obtain drops more uniform in size, add 10 drops of food of the Chemplate. This first cavity has a 1/1 concentration a.k.a. "1 part per 1".			
	nsfer 1 drop of the food coloring k to cavity 1 and then thoroughly			
4. Using the clean dropp and then rinse off the d	er, add 9 drops of clean water to ropper.	cavity 2. Use the tip of the drop	oper to gently mix the solution	
0 .	clean dropper, transfer 1 drop o	•	n cavity 2 to cavity 3.	
6. Using the clean droppe and then rinse off the d	er, add 9 drops of clean water to ropper once again.	cavity 3. Use the tip of the dropp	per to gently mix the solution	
•	of moving 1 drop of each newly conclean water followed by a gentle			
8. Construct a neat Data	Chart with the following heading	s and fill it in based on experime	ental findings.	
xample portion of data table:				
Cavity #	Concentration	Concentration	Color Observed	
1	1/1	1 part per 1	Very Dark Red	

## **Conclusions** (Answer on the Back of this Page):

- 1. At what <u>concentration</u> did you notice that the color of the dye was no longer visible? (1)
- 2. Based on your answer to guestion #1, what concentration was your detection limit for colored food dye? (1)
- 3. Which cavity has a concentration of 1 ppm? (1)
- 4. How many ppm of salt are found in a 3.5% salt solution? Show set up and work (1)
- 5. If you were the manager of a sewage treatment plant, you might be responsible for diluting treated sewage to safe levels in order to discharge them legally. If a one liter sample of the treated sewage had a concentration of 5000ppm, what would be the volume of the solution when diluted to an acceptable concentration of 5 ppm?

Formula: Concentration 1 x Volume 1 = Concentration 2 x Volume 2  $(C_1V_1 = C_2V_2)$  Show all work. (1)

**Grading:** Participation (2 points) Data Chart (3 points) Conclusion Questions (5 points)